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**Revisione:** 001 del 03/06/2015

**Composition:** BIOPLANT / BIOPLANT ELASTA bone tissue with sponge or cortical form of equine origin, composed of bone mineral and bone collagen type I. **OSTEOGEN** bone tissue with sponge form, sponge cortical tissue of equine origin, composed of bone mineral.

**Therapeutic class:** substitute of autologous or homologous bone extraction.

**Denaturation:** the denaturation process ensures the absence of rejection phenomena by the receiving body.

**Therapeutic indications:** BIOPLANT reconstruction of bone defects, even those submitted to functional charge.

**BIOPLANT ELASTA / OSTEOGEN** reconstruction of bone defects.

**Contraindications:** BIOPLANT / BIOPLANT ELASTA proved individual hypersensitivity to collagen.

**Use during pregnancy:** this product has not been tested in pregnancy conditions.

**Sterilization:** sterilized monodose product. Sterilization method by Beta ionizing radiation in 25Kgy doses.

**Shelf-life:** 5 years.

**Side effects:** as in any surgery procedure, side effects may occur.

**Caution:** This product is a mono dose, once the package is opened, any remaining cannot be reused or sterilized again. Before using it, verify the package integrity, otherwise do not use it. An inappropriate use may cause surgical nature complications.

**Preservation:** this product must be stored in a clean and dry place. The preservation temperature during the product shelflife must not be under 0°C / 32°F. DO NOT store it in refrigerated rooms or chambers.

**Transport:** product tests have proven that, during transport in normal conditions, the original packaging remains intact. Both the product and its original packaging may suffer temperature changes during short periods of time.

**Conditions of use:** pastry products packed in syringe are ready to use. Lyophilized tissue must be rehydrated before being used with physiological solution or sterilized water at ambient temperature. This product must be stabilized with osteosynthesis instruments. The product in cortical strips with code **OST- 08/09/10/15** must only be stabilized with cerclage. The products as a unit must adapt to the application area and have an intimate contact with the remaining vital bone surfaces. The product must NOT be used in an infected area. Asepsis of the product must be ensured during its manipulation.

This product MAY be used with antibiotics.

This product may NOT be used with antiinflammatories.

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